

Thank you for choosing Optilia products.

This product has been tested with Windows 7 Professional, Windows 8.1 and Windows 10. If you have problems running on other versions, please contact your dealer or send a mail to info@optilia.eu

Below you will find step by step instructions on how to install.

1. Install the driver by running "Optilia 2.0MP Camera V3.3 Setup.exe". Make sure the Flexia Definition 2.0MP or Mediscope 2.0MP camera is not connected.
2. Follow up the on-screen instructions. Restart your PC before connecting your camera.
3. Connect the camera into a free USB port of your PC. Camera will be detected automatically. Follow up the on-screen instructions if necessary.