

## 

Product No. <u>OP-006 381</u>



This fluorescence ring-light creates a flicker-free ultra violet illumination. Thanks to its high-frequency control and high quality bulb the ring-light ensures high intensity and very homogenous illumination of the objects. The UV bulb can be easily replaced by a white (daylight) bulb if required.

UV-radiation is in the range of UV-A (320 - 400 nm). Direct exposure to eyes shall therefore be avoided. UV protection glasses must be used when working with this illuminator.

## **Specification:**

UV Ring-Light	
Lamp	UV fluorescent lamp
Wave length	365 nm Max.
Dimensions	98 x 170 mm, height 37 mm
Inner Diameter	Ø 66 mm

## Including items when you order OP-006 381:

- Fluorescent UV Ring-Light Ø 66mm mounting hole
- Control unit
- Power supply and Mains cable

## **Optional accessories:**

- UV- Replacement lamp (spare part)
- White Replacement lamp



